



trails better than anything else on the market. The patent-pending technology makes blood appear to strobe like a neon sign. WATERPROOF


The IllumaTrace™ Blood Tracker is a revolutionary tool that helps hunters detect and identify blood















**CHARGING INSTRUCTIONS** 

Charging time may take up to 3 hours.

**Solid Red:** Charging

Solid Green: Fully Charged

**Steady Orange:** Battery more than 25% Pulsing Orange: Battery less than 25%

**Muma**Trace<sup>®</sup>

NEBO


source with a USB-A port.

WHEN CHARGING

WHEN IN USE

Remove the battery from the flashlight and insert the USB-C end of the included charging cable into

the USB-C port on the battery. Next, insert the USB-A end of the charging cable into a power



**OPERATION** 

Traditional flashlight. IllumaTrace™ Modes(S1, S2, S3)

**MODES** 


White Mode(W)

strobe rate that works best for your environment.

WHITE MODE **ILLUMATRACE™** 

**LUMENS** 

Choose between S1, S2, and S3 to determine the





2h 30m 262 ft

**BEAM DISTANCE** 

246 ft

**RUN TIME** 

3h

If you experience a problem, please do not WARRANT return to your local store. Contact our Customer Service team at 800-255-6061 or **NEBOLIGHTS.com** 

NEBOLIGHTS.com

For more information, detailed video

instructions, and warranty registration, visit