



trails better than anything else on the market. The patented technology makes blood appear to strobe like a neon sign.

The IllumaTrace® Blood Tracker is a revolutionary tool that helps hunters detect and identify blood
















WHEN CHARGING A small light near the positive terminal of the

Solid Green: Fully Charged

Solid Red: Charging

BATTERY LEVEL INDICATOR

battery serves as the charging status indicator.

(Charing time may take up to 3 hours.)

The power button light serves as the

Remove the battery from the flashlight and insert the USB-C end of the included charging cable into

the USB-C port on the battery. Next, insert the USB-A end of the charging cable into a power

CHARGING INSTRUCTIONS

source with a USB-A port.

battery level indicator while the product is in use (any mode). **Steady Orange:** Battery more than 25% Flashing Orange: Battery less than 25%

NumaTrace



once. It will illuminate and then automatically turn

White Mode (W)

off after 5 seconds.

Choose between S1, S2, and S3 to determine the strobe rate that works best for your environment.

turned on accidentally or while in storage, thereby

checked in this mode by pressing the power button

conserving battery life. The battery level can be

flashlight mode. IllumaTrace® Modes(S1, S2, S3)

This mode is a non-strobing, traditional

WHITE MODE **ILLUMATRACE®**



170



WHITE

If you experience a problem, please do not return to your local store. Contact our Customer Service team at 800-255-6061 or NEBOLIGHTS.com

instructions, and warranty registration, visit

For more information, detailed video

NEBOLIGHTS.com

2h 30m

262 ft